

Comparison of Automated Technology (SQA-V) and Manual Technology (Microscope) for Routine Semen Analysis

Sperm Concentration	Microscope per WHO		Automated SQA-V	
Sample preparation & counting methodology	Mix the sample. Use a hemacytometer or a Makler counting chamber. Follow counting chamber instructions. Count 200 sperm cells twice, verify that the difference between the two counts is within the acceptable statistical range. If not, the sample must be rerun.		Mix the sample. Fill a capillary and follow the on-screen instructions for automated testing .	
Technology	Manual count under the microscope. Accuracy and precision are dependent on observer proficiency and experience.		Automated System based on the objective electro-optical method (spectrophotometry).	
Source of errors	Subjective nature of manual counting leads to high inter- and intra-operator discrepancies. Test results will vary depending on the type of counting chamber used.			
Sample size tested	10 µl	The overall specimen is poorly represented because of the small sample size.	650 µl	The overall specimen is adequately represented because of the large sample size.
Sample fraction tested (2-5 ml sample)	0.2% - 0.5%		13% - 32.5%	
Number of cells measured	200 x 2		Millions of cells	
Results	Counted manually - based on the type of counting chamber used. Results include counting and subjective errors.		Fully objective automated test results archived in a database.	
Average testing time	15-20 min.		75 seconds for all semen parameters.	
Accuracy	The actual coefficient of variation > 9-20% . The acceptable coefficient of variation is up to 10% if counting 200 sperm cells (WHO'99 manual). The difference between duplicate counts of 200 sperm cells can reach 14%.		The actual coefficient of variation < 2% .	

Sperm Motility	Microscope per WHO		Automated SQA-V	
Sample preparation & counting methodology	Mix the sample. Use a standard slide or a Makler counting chamber. Follow counting instructions of the WHO and differentiate cells according to a velocity grading system. Count 200 sperm cells twice, verify that the difference between the two counts is within the acceptable statistical range. If not, the sample must be rerun.		Mix the sample. Fill a capillary and follow the on-screen instructions for automated testing	
Technology	Manual count under the microscope. Accuracy and precision are dependent on observer proficiency and experience.		Automated System based on the objective electro-optical method (signal processing of the light modulations caused by the motile sperm cells).	
Source of errors	Subjective nature of manual counting leads to high inter- and intra-operator discrepancies. Motility is frequently over-estimated.			
Sample size tested	10 µl		50 µl	
Sample fraction tested (2-5 ml sample)	0.2%-0.5%	The overall specimen is poorly represented because of the small sample size.	1%-2.5%	The overall specimen is adequately represented because of the large sample size.
Number of cells measured	200 x 2		Thousands	
Results	Manual count and calculation based on the number of cells for each type of velocity grading system: a, b, c, d.		Fully objective automated test results archived in a database.	
Average testing time	15-20 min.		75 seconds for all semen parameters.	
Accuracy	The actual coefficient of variation > 9-20% . The acceptable coefficient of variation is 14% - 50% if counting 200 sperm cells (WHO'99 manual). The percentage difference between duplicate counts of 200 sperm cells can reach 10 % (absolute units).		The actual coefficient of variation < 3% .	

Sperm Morphology	Microscope per WHO		Automated SQA-V	
Sample preparation & counting methodology	Mix the sample. Follow WHO fixing and staining instructions for Papanicolaou method. Differentiate morphological abnormalities using WHO or Kruger strict criteria. Count 200 sperm cells twice, verify that the difference between the two counts is within the acceptable statistical range. If not, the sample must be rerun.		Mix the sample. Fill a capillary and follow the on-screen instructions for automated testing	
Technology Source of errors	Manual assessment and count under the microscope. Accuracy is highly dependent on the proficiency and experience of the technician/MD. High inter- and intra-operator variability due to the subjective nature of the test; requiring the ability to accurately assess and differentiate sperm cells with morphological abnormalities.		Automated test. Calculated parameter (proprietary algorithm) based on the correlation between progressive motility and normal morphology.	
Sample size tested Number of cells measured	10 µl 200 x 2	The overall specimen is poorly represented because of the small sample size.	50 µl Thousands	The overall specimen is adequately represented because of the large sample size.
Results Average testing time	Manual grading system based on the assessment and classification of individual spermatozoa characteristics and abnormalities (head, neck, mid-piece, etc according to WHO). 15-20 min.		Fully automated test results archived in a database. 75 seconds for all semen parameters.	
Accuracy	The actual coefficient of variation > 9-40% . The acceptable coefficient of variation is 14% - 50% if counting 200 cells (WHO'99 manual). The percentage difference between duplicate counts of 200 sperm cells can reach 10 % (absolute units)		The actual coefficient of variation < 3% .	