



LABMATE Soft
single channel pipette



SCIENCE AT YOUR FINGERTIPS

LABMATE Soft

single channel pipette

For years LABMATE variable volume pipettes have been the highlight of the HTL product range. Their outstandingly accurate and precise pipetting mechanism, well balanced and robust handle made LABMATE the pipette of choice for many laboratory professionals worldwide. The most recent version of this popular pipette - **LABMATE Soft** - combines so appreciated traditional design with the innovative thumb-friendly spring system and new volume setting pushbutton, making the daily pipetting routine an effortless and reliable experience.



RECALIBRATION SYSTEM
allows for quick but accurate pipette recalibration



VOLUME SETTING
using the pushbutton or thumbwheel



FINGER - HOOK
takes the weight, for a more relaxed grip

REDESIGNED SOFT SPRING SYSTEM
requires much lower pressures for pipetting



LARGER VOLUME DISPLAY
with easy-to-read numbers

UV RESISTANT BODY
can be safely sterilized with UV radiation - even prolonged exposure



HEIGHT ADJUSTABLE EJECTOR
to accommodate virtually all brand of tips

HIGHLY DURABLE SHAFT
constructed from tough PVDF guarantees highest mechanical and chemical resistance



Dual volume adjustment system

In order to meet preferences of all pipette users **LABMATE Soft** pipettes are equipped with a dual adjustment system. Pipette volume can be set precisely with one hand by the turn of the volume setting knob, whilst for rewinding of the volume counter the pipetting pushbutton can be applied. The top part of the pushbutton is rotational to prevent accidental changing of the pipette volume.



Pipette stands

In order to ensure convenient storage of pipettes, HTL offers a wide range of single and multiple stands.



5446

CAROUSEL



5484

ROTATING



5447

MULTIPLE



5480

LINEAR

HTL pipette tips

HTL Pipette Tips are designed to fit a wide variety of single and multichannel pipettes. Flexible walls and the diameter of tips collars ensure a secure fit with less force required to load and eject. All pre-sterilized tips are exposed to e-beam radiation for complete security without degrading the plastic or performance. Tips are manufactured in one of the most technically advanced plastic disposable production facility in clean room conditions. HTL Pipette Tips are certified to be free from bio active contaminants as RNase, DNase or pyrogens. HTL tips are offered in standard, sterile and filtered versions in all popular packaging formats:



BULK



RACK



STACK
RACK



REFILL
SYSTEM

Distributed by:

All rights reserved. Product described in this brochure is subject to availability and technical modification. Errors excepted. PZ HTL S.A. reserves the right to improve, enhance or otherwise modify its products without prior notification.

