



# **SWIFTPET PRO**

**PIPETTE CONTROLLER**



***suitable for 0.5 ml to 100 ml pipettes***

***protected by filter***

***ergonomically shaped handle***

***LCD display indicates the level  
of battery charge***



SWIFTPET PRO - lightweight and ergonomically shaped motorized pipette controller designed for work with glass or plastic serological pipettes in the 0.5 - 100 ml range. Easily accessible switches allow choosing different operation modes depending on the volume of a pipette and liquid viscosity. Switches enable selection of HIGH/LOW speed and also gravity/blow-out mode dispense. Additionally aspirating and dispensing speed is controlled through two comfortable finger triggers. When fully charged, the NiMH batteries enables up to eight hours of continuous use. LCD display indicates setting modes and level of battery charge. In order to protect the device against overfilling the pipette controller is equipped with a 0.2 µm PTFE filter. The filter, the pipette holder and the nosepiece are easily removed for autoclaving. The SWIFTPET PRO is also supplied with the universal battery charger, an additional hydrophobic 0.2 µm filter and the bench/wall mounting charging stand.



## Ordering information



Model	Cat. no.
SWIFTPET PRO with charging stand	0390

Charging stand	29056
Pipette holder	29054
Nosepiece	29053
Micro filter 0.20 µm	9023
Micro filter 0.45 µm	9024
Universal charger 9V type: EU, US, UK, AU	9010

Distributed by:

PZ HTL S.A. ■ DANISZEWSKA 4 ■ 03-230 WARSAW ■ POLAND  
 Tel.: +48 22 492 19 00 ■ Fax: +48 22 492 19 93 ■ E-mail: htl.sales@corning.com  
 www.htl.com.pl

